

New Generation Access System

- ▶ Features superior fixation
- ▶ Compatible with Kii obturators
- ▶ Advanced seal design

Kii[®] Advanced Fixation

Access System



Kii® Advanced Fixation Systems



Retention Disk

*Allows fixation of the trocar at specific depths
Color-coded to indicate cannula size*

Fixation Balloon

*Maintains superior cannula retention without barbs
Eliminates unintentional displacement of trocar
Maximizes working space*

Available with Fios® first entry, shielded bladed and optical obturators

Advanced Seal Design

*Offers high flow insufflation and exceptional durability
Facilitates simple seal detachment for rapid desufflation
Ensures smooth instrument exchanges with minimal drag*



SIZE & LENGTH	REORDER NO.			
	FIOS FIRST ENTRY	OPTICAL ACCESS SYSTEM	SHIELDED BLADED ACCESS SYSTEM	CANNULA & SEAL
5 x 75mm	CFF05	N/A	N/A	CFS03
5 x 100mm	CFF03	CFR03	CFB03	CFS02
5 x 150mm	CFF01	CFR01	CFB01	CFS01
11 x 100mm	CFF33	CFR33	CFB33	CFS12
11 x 150mm	CFF31	CFR31	CFB31	CFS11
12 x 100mm	CFF73	CFR73	CFB73	CFS22
12 x 150mm	CFF71	CFR71	CFB71	CFS21

All models consist of six units per box, except for the cannula and seal models which are twelve per box. Please contact your Sales Representative for additional configurations.

Applied Medical
A New Generation Medical Device Company

Rancho Santa Margarita, CA 92688
Tel 949.713.8300
Fax 949.713.8400
www.appliedmedical.com

To place an order, or for more information, call 800.282.2212



© 2011 Applied Medical Resources Corporation. All rights reserved.
Applied Medical, the Applied Medical logo design, ProTerra, the ProTerra logo design and marks designated with a ® are trademarks of Applied Medical Resources Corporation, registered in one or more of the following countries: Australia, Canada, Japan, the United States and/or the European Union. SC01399E

