


Advanced Transanal Access

- ▶ Facilitates minimally invasive transanal surgery including TEM/TAMIS
- ▶ Saves valuable O.R. time with simple set-up
- ▶ Provides continuous access and visualization

GeIPOINT® Path

Transanal Access Platform






GelPOINT® Path

GelSeal® Cap

- Provides a flexible fulcrum for unmatched triangulation of standard laparoscopic instrumentation
- Maintains pneumorectum for continuous access
- Incorporates smoke evacuation capabilities for continuous visualization
- Facilitates large specimen removal with simple detachment from access channel



Actual size

Access Channel

- Saves valuable OR time with simple placement and set-up
- Offers 4cm atraumatic retraction for enhanced exposure and access
- Includes retention flange and suture ties for maximum security throughout procedure

Self-Retaining Ports

- Accommodates 5-10mm standard instrumentation, 5mm angled instrumentation and curved suture needles
- Offers greater freedom of movement and less instrument collision due to low-profile design



REORDER NO.	DESCRIPTION	ACCESS CHANNEL DIAMETER/LENGTH	QUANTITY
CNO11	GelPOINT Path	4cm x 4.5cm	1/BOX

GELPOINT PATH COMPONENTS

	SIZE	QUANTITY
GelSeal Cap with Insufflation Port		1
Access Channel with Introducer	4cm x 4.5cm	1
Self-Retaining Ports	5-10mm	3
Obturator	5mm	1



Rancho Santa Margarita, CA 92688
 Tel 949.713.8300
 Fax 949.713.8400
www.appliedmedical.com
www.gelpoint.com

To place an order, or for more information, call 800.282.2212



© 2011 Applied Medical Resources Corporation. All rights reserved.
 Applied Medical, the Applied Medical logo design, ProTerra, the ProTerra logo design and marks designated with a ® are trademarks of Applied Medical Resources Corporation, registered in one or more of the following countries: Australia, Canada, Japan, the United States and/or the European Union. SC01482C

