


EVACUATE SMOKE AT THE SOURCE WITH **EPIX[®] ELECTROSURGICAL PROBES**

The Epix probes deliver monopolar electrosurgery while effectively removing smoke at the source. The gradual on-demand button allows for smoke evacuation while minimizing loss of insufflation, which is ideal when operating at lower pressures.



EPIX ELECTROSURGICAL PROBES



Epix Electrosurgical Probes with Smoke Evacuation, Angled

Reorder No.	Tip Design	Size	Qty/Box
CW001*	Spatula	5mm x 42cm	1
CW002*	L-hook	5mm x 42cm	1
CW003*	Needle	5mm x 42cm	1



Epix Electrosurgical Probes with Smoke Evacuation, Straight

Reorder No.	Tip Design	Size	Qty/Box
CW004**	Spatula	5mm x 42cm	1
CW005**	L-hook	5mm x 42cm	1
CW006**	Needle	5mm x 42cm	1

* Models CW001, CW002 and CW003 are compatible with Applied Medical GelPOINT® sleeves and Applied Medical Kii® 12 x 60-100mm and 15 x 100-150mm trocars.

** Models CW004, CW005 and CW006 are compatible with Applied Medical GelPOINT sleeves and Applied Medical Kii 5mm or larger trocars (if using non-Kii trocars, a larger port size will be needed).

Visit www.appliedmedical.com/epix
for more information

Devices listed may not be approved in all markets.
Please contact your Field Implementation team member for more information on availability.

© 2020 Applied Medical Resources Corporation. All rights reserved.
Applied Medical, the Applied Medical logo design and marks designated with a ® are trademarks of Applied Medical Resources Corporation, registered in one or more of the following countries: Australia, Canada, Japan, South Korea, the United States, and/or the European Union. All of the other trademarks used herein are the property of their respective owners. 300841-EN-USA-A