



Epix[®]

Laparoscopic Instrumentation

Epix

Universal Clip Applier



Advanced Clip Delivery System

Offers reliable, smooth device actuation.

34cm Shaft

Maximizes instrument reach.

Ergonomic Handle Design

Provides optimal control during clip loading and closure.



Low-Profile Jaw

Allows for maximum visualization and clearance around the tissue to be ligated.


Allows for greater clip retention through the application of a true M/L clip.



True M/L clip

Superior Clip Retention

Features TiGold® clip technology, ensuring reliable closure and secure retention.



APPROXIMATELY
100%
GREATER RETENTION
COMPARED TO 5mm CLIP APPLIERS

CLIP RETENTION FORCE (pounds)

Note: Data is on file.

Epix

Electrosurgical Probes with Smoke Evacuation

Probe Tips

Provide precise delivery of monopolar energy.




CW001



CW002



CW003




Angled Shaft

Improves instrument reach and visibility
Reduces instrument clashing

Built-In Smoke Evacuation

Evacuates smoke at the push of a button while minimizing loss of insufflation.



Sliding Sheath

Advances for protection during instrument insertion.
Retracts to expose electrosurgical tip.

Epix

Scissors, Graspers, and Dissectors

Epix Scissors

High-performance blades sustain sharpness throughout the duration of the procedure. Low-profile design eliminates protrusions, enabling full insulation down to the blades.

Epix Graspers

Clinically proven *latis*® pads gently engage tissue, reducing the pressure necessary to maintain traction.

Handle design incorporates the ability to grasp and lock tissue as desired, ensuring precise control.

Epix Dissectors

Curved jaws enhance the grasping, mobilizing and cauterizing of tissue.

Tapered, serrated jaws provide traction to engage tissue.



Epix Scissors

Reorder No.	Size	Qty/Box
<i>Laparoscopic</i>		
CB030	5mm x 35cm	10
CB040	5mm x 45cm	10



Epix *latis* Grasper

Reorder No.	Size	Qty/Box
<i>Laparoscopic</i>		
C4130	5mm x 35cm	10
C4140	5mm x 45cm	10



Epix Babcock Grasper

Reorder No.	Size	Qty/Box
<i>Laparoscopic</i>		
C4230	5mm x 35cm	10
C4240	5mm x 45cm	10



Epix Dissectors

Reorder No.	Size	Qty/Box
<i>Laparoscopic</i>		
CY010	5mm x 35cm	10
CY020	5mm x 45cm	10



Epix Universal Clip Applier

Reorder No.	Size	Qty/Box
-------------	------	---------

(Contains 20 M/L TiGold clips)

CA500*	5mm x 34cm	3
--------	------------	---



Epix Electrosurgical Probes with Smoke Evacuation, Angled

Reorder No.	Tip Design	Size	Qty/Box
-------------	------------	------	---------

CW001**	Spatula	5mm x 42cm	1
---------	---------	------------	---

CW002**	L-hook	5mm x 42cm	1
---------	--------	------------	---

CW003**	Needle	5mm x 42cm	1
---------	--------	------------	---

* Model CA500 is compatible with Applied Medical Kii® 5mm and larger trocars.

** Models CW001, CW002 and CW003 are compatible with Applied Medical GelPOINT sleeves and Applied Medical Kii 12 x 60-100mm and 15 x 100-150mm trocars.

Visit www.appliedmedical.com/epix for more information.

Devices listed may not be approved in all markets.

Please contact your Field Implementation team member for more information on availability.

© 2019 Applied Medical Resources Corporation. All rights reserved.

Applied Medical, the Applied Medical logo design and marks designated with a ® are trademarks of Applied Medical Resources Corporation, registered in one or more of the following countries: Australia, Canada, Japan, South Korea, the United States, and/or the European Union. All of the other trademarks used herein are the property of their respective owners. 143331-EN-USA-C

