



Epix[®]
Universal Clip Applier

ONE SOLUTION FOR LIGATION



Low-Profile Jaw

Allows for maximum visualization and clearance around the tissue to be ligated.

Delivers a true M/L clip through a 5mm shaft when the device is used in tandem with a Kii 5mm trocar.



34cm Shaft

Maximizes instrument reach.

Advanced Clip Delivery System

Offers reliable, smooth device actuation.

Ergonomic Handle Design


Provides optimal control during clip loading and closure.



TiGold® Clip Technology

Ensures reliable closure and secure clip retention.

True M/L clip



CLIP RETENTION FORCE (pounds)

Note: Data is on file.



Epix Universal Clip Applier

Reorder No.

Qty/Box

Epix universal clip applier (contains 20 M/L TiGold clips)

CA500*

3

* The Epix universal clip applier is compatible with Applied Medical Kii® 5mm trocars and larger trocars.

Visit [appliedmedical.com/epix](https://www.appliedmedical.com/epix) for more information.

Please contact your Applied Medical representative for more information on availability. This information is intended for dissemination exclusively to healthcare professionals and is not intended to replace labeling and Instructions for Use (IFU). Please refer to the IFU for the indications, contraindications, warnings, precautions, instructions and other information.

© 2021 Applied Medical Resources Corporation. All rights reserved. Applied Medical, the Applied Medical logo design, and marks designated with the symbol © are trademarks of Applied Medical Resources Corporation, registered in one or more of the following countries: Australia, Canada, Japan, South Korea, the United States, the United Kingdom, and/or the European Union.

11/2021 194619-EN-USA-B

