



Epix[®] Universal

Clip Applier

ONE SOLUTION FOR LIGATION



Low Profile Jaw

Allows for maximum visualization and clearance around the tissue to be ligated

Allows for greater clip retention through the application of a true M/L clip



34cm Shaft

Maximizes instrument reach

Advanced Clip Delivery System

Offers reliable, smooth device actuation

Ergonomic Handle Design


Provides optimal control during clip loading and closure



True M/L clip

Superior Clip Retention

Features TiGold® clip technology, ensuring reliable closure and secure retention





Epix Universal Clip Appliers

Reorder No.

Qty/Box

Epix Universal (contains 20 M/L TiGold clips)

CA500*

3/Box

* Epix Universal clip applier is compatible with Applied Medical Kii® 5mm and larger trocars.

Visit www.appliedmedical.com/epix
for more information

© 2017 Applied Medical Resources Corporation. All rights reserved.
Applied Medical, the Applied Medical logo design and marks designated with a © are trademarks of Applied Medical Resources Corporation, registered in one or more of the following countries: Australia, Canada, Japan, South Korea, the United States and/or the European Union. All of the other trademarks used herein are the property of their respective owners. SC02497A

