



GelPOINT[®]

V-Path

Transvaginal Access Platform

GelPOINT V-Path

Transvaginal Access Platform

GelSeal® Cap


Provides a flexible fulcrum for unmatched triangulation of standard laparoscopic instrumentation

Maintains insufflation for continuous access

Includes two valves to facilitate smoke evacuation and accommodate standard laparoscopic insufflation tubing

Facilitates large specimen extraction with simple detachment from Alexis® retractor

*Two GelSeal Cap
sizes available*



GelPOINT V-Path Transvaginal Access Platform


Enables access for minimally invasive instruments to perform transvaginal natural orifice transluminal endoscopic surgery (vNOTES) procedures



Self-Retaining Sleeves

Accommodate 5-12mm standard laparoscopic instrumentation

Offer greater freedom of movement and less instrument collision due to low-profile design



Alexis Retractor

Accommodates varying incision sizes

Offers 360° of atraumatic retraction for enhanced exposure and access

Three sizes available to meet procedural needs



GelPOINT V-Path Transvaginal Access Platform

Reorder No.	Alexis Retractor Inner Ring Size	Qty/Box
C2A11	11cm	1/Box

Components

(1) GelSeal Cap
(1) Alexis Retractor
(3) 10mm Sleeves
(1) 12mm Sleeve
(1) Obturator
(1) Instrument Shield



Reorder No.	Alexis Retractor Inner Ring Size	Qty/Box
C2A12	9.5cm	1/Box

Components

(1) GelSeal Cap
(1) Alexis Retractor
(3) 10mm Sleeves
(1) 12mm Sleeve
(1) Obturator
(1) Instrument Shield



Reorder No.	Alexis Retractor Inner Ring Size	Qty/Box
C2A15	7cm	1/Box

Components

(1) GelSeal Cap
(1) Alexis Retractor
(3) 10mm Sleeves
(1) 12mm Sleeve
(1) Obturator

Visit www.appliedmedical.com
for more information

Devices listed may not be approved in all markets.

Please contact your field implementation team member for more information on availability.

© 2019 Applied Medical Resources Corporation. All rights reserved.

Applied Medical, the Applied Medical logo design and marks designated with a ® are trademarks of Applied Medical Resources Corporation, registered in one or more of the following countries: Australia, Canada, Japan, South Korea, the United States, and/or the European Union.

140625-EN-EUR-A

