






Kii[®]
Abdominal Access Systems

Kii Fios[®]

First Entry Access System

THE ONLY INSUFFLATING OPTICAL PORT



| Kii Fios Access System | Insufflation Needle | 5mm Access System | 12mm Access System |
|--|---|---|---|
| 3mm |  |  |  |
|  | | | |
|  | | | |
| | | | Minimal Fascial Defect |
| | | | True nonbladed access allows for reapproximation of tissue layers with minimal fascial defect once the trocar is removed. |
| | | | Minimal Penetration with Rapid Insufflation |
| | | | Only 3mm of trocar penetration is needed to initiate insufflation. |

Kii

Balloon Blunt-Tip Access System

INNOVATIVE BLUNT-TIP TECHNOLOGY



Suture Ties

Maximize instrument reach and facilitate port-site closure.

Non-Latex Balloon

Provides stable anchoring throughout lengthy procedures while maximizing working space and visibility.



Now in the 5mm size — the smallest Hasson trocar available!

Gel Cone

Enhances instrument articulation.


Optical Blunt-Tip Obturator

Facilitates direct visualization during trocar placement.

Kii

Advanced Fixation Technology

SUPERIOR RETENTION




Retention Disk

Reduces migration of the trocar sleeve into the abdominal cavity.



Non-Latex Balloon

Provides unmatched retention.
Maximizes working space and visibility.
Provides pressure on the incision site, potentially reducing port-site bleeding.
Maintains stability throughout procedures.



Kii

Low-Profile Access System

SLIM, LIGHTWEIGHT DESIGN



Low-Profile Seal Housing

Optimizes external working space, minimizing interference with adjacent ports.



Low-Profile Kii Seal



Standard 5mm Kii Seal

Z-Thread Sleeve with Extended Threads

Maintains retention at shorter depths, maximizing internal working space and visibility.

Kii Low-Profile Access System

| Reorder No. | | | Size |
|----------------|---------------|---------------|----------------------|
| <i>Optical</i> | <i>Bladed</i> | <i>Sleeve</i> | |
| CTR21 | CTB21 | CTS28 | 5 x 55mm |
| CTR23 | CTB23 | N/A | 5 x 55mm, dual pack |
| CTR22 | CTB22 | CTS29 | 5 x 100mm |
| CTR24 | CTB24 | N/A | 5 x 100mm, dual pack |


Kii

Shielded Bladed Access System

SMOOTH, CONTROLLED ENTRY

Intuitive Switch Design

Provides a clear indication of shield position.



Deactivated



Activated

Kii Shielded Bladed Access System

| Reorder No. | | Size |
|-------------|--------------------------|-----------------|
| | <i>Advanced Fixation</i> | <i>Z-Thread</i> |
| N/A | CTB05 | 5 x 55mm |
| CFB03 | CTB03 | 5 x 100mm |
| CFB01 | CTB01 | 5 x 150mm |
| N/A | COQ61 | 8 x 100mm |
| CFB33 | CTB33 | 11 x 100mm |
| CFB73 | CTB73 | 12 x 100mm |
| CFB71 | CTB71 | 12 x 150mm |

Kii

Optical Access System

UNPARALLELED OPTICAL CLARITY



Tissue Reapproximation



Initial
Withdrawal



Midway



Complete
Withdrawal

Kii Optical Access System

| Reorder No. | | Size |
|--------------------------|-----------------|---------------------|
| <i>Advanced Fixation</i> | <i>Z-Thread</i> | |
| N/A | CTR05 | 5 x 55mm |
| CFR03 | CTR03 | 5 x 100mm |
| CFR01 | CTR01 | 5 x 150mm |
| N/A | COQ19 | 8 x 100mm |
| CFR33 | CTR33 | 11 x 100mm |
| CFR34 | CTR34 | 11 x 100mm, handled |
| N/A | COR83 | 12 x 60mm |
| CFR73 | CTR73 | 12 x 100mm |
| CFR74 | CTR74 | 12 x 100mm, handled |
| CFR71 | CTR71 | 12 x 150mm |
| N/A | CTR72 | 12 x 150mm, handled |
| N/A | COR37 | 15 x 100mm |
| CFR39 | COR39 | 15 x 150mm |



Kii Fios First Entry Access System

Reorder No. Size

Advanced Fixation Z-Thread

| | | |
|-------|-------|---------------------|
| CFF05 | N/A | 5 x 75mm |
| CFF03 | CTF03 | 5 x 100mm |
| CFF04 | CTF04 | 5 x 100mm, handled |
| CFF01 | CTF01 | 5 x 150mm |
| CFF33 | CTF33 | 11 x 100mm |
| CFF34 | CTF34 | 11 x 100mm, handled |
| CFF75 | N/A | 12 x 75mm |
| CFF73 | CTF73 | 12 x 100mm |
| CFF74 | CTF74 | 12 x 100mm, handled |
| CFF71 | CTF71 | 12 x 150mm |



Kii Balloon Blunt-Tip Access System

Reorder No. Size

Blunt-Tip

| | |
|-------|------------|
| COQ20 | 5 x 100mm |
| COR47 | 12 x 100mm |
| COR50 | 12 x 130mm |



Kii Sleeve with Advanced Fixation

Reorder No. Size

Advanced Fixation

| | |
|-------|------------|
| CFS03 | 5 x 75mm |
| CFS02 | 5 x 100mm |
| CFS01 | 5 x 150mm |
| CFS12 | 11 x 100mm |
| CFS22 | 12 x 100mm |
| CFS21 | 12 x 150mm |

Visit www.appliedmedical.com/trocars
for more information.