



Insufflation Needle

A large teal graphic element consisting of a diagonal line that starts from the top left and extends towards the bottom right, creating a triangular shape on the left side of the page.

ENABLES HIGH FLOW INSUFFLATION

13 gauge percutaneous needle

Tactile and visual indicators

Ergonomic housing design

Insufflation Needle



Spring Loaded Blunt Stylet with Red Indicator Button

Provides tactile and visual feedback when stylet is retracted

13 Gauge High Flow Needle

Offers faster insufflation than traditional needles, reducing OR time



Ergonomic Housing Design

Ensures control during needle insertion





Insufflation Needle

Recorder No.	Size	Quantity
C2201	13 Gauge, 120mm	20/Box
C2202	13 Gauge, 150mm	20/Box

Visit www.appliedmedical.com/laparoscopic
for more information

© 2018 Applied Medical Resources Corporation. All rights reserved.
Applied Medical, the Applied Medical logo design and marks designated with a ® are
trademarks of Applied Medical Resources Corporation, registered in one or more of the
following countries: Australia, Canada, Japan, the United States, and/or the European Union.
SC01450D

