



VOYANT[®]
INTELLIGENT ENERGY SYSTEM

Voyant

Intelligent Energy System



Rather than relying on only one variable, the Voyant system intelligently measures the tissue and its changing properties and applies the optimal amount of energy to create a permanent, fused seal.

Intelligence

The Voyant system introduces a unique advancement by embedding the intelligence in each device. This allows an unlimited potential for future instruments that provide precise solutions for specific clinical needs.




Voyant 5mm Fusion



Jaw Features

- 20mm seal length
- 18mm cut length
- Straight, blunt jaw



Voyant Open Fusion



Jaw Features

- 40mm seal length
- 38mm cut length
- Curved, blunt jaw



Voyant Fine Fusion



Jaw Features

- 17mm seal length
- 15mm cut length
- Curved jaw with dissecting tips





Laparoscopic Devices

Reorder No.	Description	Size	Qty/Box
EB010	Voyant 5mm Fusion	5mm x 37cm	6/Box
EB011	Voyant 5mm Fusion	5mm x 44cm	6/Box



Open Devices

Reorder No.	Description	Size	Qty/Box
EB030	Voyant Fine Fusion	13.5cm	6/Box
EB040	Voyant Open Fusion	20cm	6/Box



Generator

Reorder No.	Description	Qty/Box
EA010	Voyant Electrosurgical Generator	1/Box



Accessories

Reorder No.	Description	Qty/Box
EX050	Voyant Cart	1/Box
EX051	Voyant Cart Cover	1/Box

Visit www.appliedmedical.com/Voyant
for more information